

Upgraded Model

LMSPV325



Visible Spectrophotometer

Salient Features

✦ MULTIPLE RESULTS READOUT

Can display wavelength, absorption and transmittance with 5 results per screen. It also has a memory store of up to 200 results.

✦ AUTO SETTING WAVELENGTH

Users set wavelength automatically through arrow keys to avoid operation errors.

✦ STANDARD CURVE

Can set up various standard curves according to customer's solutions and find the concentration of unknown solutions.

✦ DATA OUTPUT

LMSPV325 is equipped with USB port to connect with a PC to display spectrum scanning, kinetics and multi wavelength testing results on the screen.

✦ 2.5" LCD SCREEN

Equipped with a 2.5 LCD screen to give a clear display of standard curves and groups of results.

Visible Spectrophotometer

HSN Code : 90273020

SPECIFICATIONS	Model	LMSPV325
	Make	LABMAN
	Optical System	Single beam, Grating 1200 lines/mm
	Wavelength Range	325-1000nm
	Bandwidth	4nm
	Wavelength Accuracy	±1nm
	Wavelength Repeatability	±0.5nm
	Wavelength Setting	Auto
	Photometric Accuracy	±0.5%T
	Photometric Repeatability	±0.3%T
	Photometric Range	-0.3-3A,0-200%T
	Stray Light	≤0.3%T
	Stability	± 0.002A/h @ 500nm
	Display	LCD
	Detector	Silicon Photodiode
	Standard Cell Holder	4-position 10mm cell holder
	Light Source	Tungsten Lamp
	Output	USB Port & Parallel Port(Printer)
Power Supply	AC 220V/50Hz	
Dimension (LxWxH) mm	420x280x180	
Weight	10kg	

★ Specifications are subject to change without notice.

STANDARD ACCESSORIES	Software	1 pc
	10mm Glass Cuvette	4 pcs.
	Power Cord	1 pc
	Instruction Manual	1 pc.
	Black Block	1 pc.
	Dust Cover	1 pc.

