

MSZ

Stereo Zoom Microscope

Biological & Industrial Use

- Zoom ratio 1 : 7
- Parfocality - one time focus for all zoom magnifications
- Ergonomic design for ease of operation
- Both transmitted and reflected illumination option
- 6v 15w halogen top lamp / LED
- 5w fluorescent bottom lamp / LED



Model	Head	Zoom Ratio	Objective Zoom Range	Eyepiece	Working Distance	Illumination
MSZ Bi	60° Binocular	1 : 7	0.65x - 4.5x	WF 10x (FN. 22)	100 mm	6v 15w Halogen / 5w Fluorescent / LED
MSZ Tr	60° Trinocular	1 : 7	0.65x - 4.5x	WF 10x (FN. 22)	100 mm	6v 15w Halogen / 5w Fluorescent / LED

Magnus