



Salient Features

- ✦ Based on E.M.F.C. Technology
- ✦ A built in under hook for below weighing
- ✦ RS232C Interface
- ✦ Tare range upto max. Capacity of the Balance
- ✦ Automatic External Calibration
- ✦ Multiple Weighing Units
- ✦ S.S. Weighing Pan, Level Indicator, Adjustable Fit
- ✦ Piece Counting

High Precision Balance

HPB3000

HSN Code - 90160020

MODEL	HPB3000
Make	WENSAR
Capacity	3000g
Readability	0.01g (10mg)
Repeatability (±)	0.01g
Linearity (±)	0.03g
Display	VFD (Vacuum Fluorescent Display)
Pan Size	150 mm Ø
Response Time	3 - 4 Sec
Calibration	Automatic External
Tare Range	Full
Operating Temperature	10°C to 30°C
Interface	RS232C
Power Supply	AC 220V / 50-60 Hz
Dimension (LxWxH)	530 x 435 x 345 mm
Weight (app.)	6.9kg

SPECIFICATIONS

★ Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, SS PAN, PAN Supporter, Adaptor, Cloth, Instruction Manual, Test Certificate.