

# DIGITAL LEAK TEST APPARATUS

75



NEW

## TECHNICAL SPECIFICATIONS

- |                             |  |
|-----------------------------|--|
| ✓ Type                      | GMP Model, In-built Test Unit                        |
| ✓ Time Control & Range      | Digital, Programmable, 00.01 sec to 99min59sec       |
| ✓ Test Vacuum Range         | Adjustable up-to 650 mmHg                            |
| ✓ Vacuum Display            | 2.5" SS Vacuum Gauge                                 |
| ✓ Vacuum Control & Release  | By Single Front Panel Knob                           |
| ✓ Audio Alarm at completion | Yes  |
| ✓ Vacuum Pump               | Oil free, MV 27" +- 5%, 1/4HP Motor                  |
| ✓ Test Unit Size            | 200mm, 250mm & 300mm (any one)<br>(Specify in Order) |

## SALIENT FEATURES

- ✓ Programmable Time With Digital Display
- ✓ Auto Stop After Completion Of Set Time With Alarm
- ✓ Precise Control of Vacuum
- ✓ Heavy Duty Oil-free Vacuum Pump
- ✓ 1mm Thick SS-304 cabinet

