

Since 1954  
**WESWOX**<sup>®</sup>  
OPTIK

# ADVANCE ROTARY MANUAL MICROTOME

Introducing  
**World's First Automated Microtome**  
with Feather Touch LCD Screen Panel.



## ARMT-1090A

**WESWOX** Advance Manual Microtome is Ergonomically designed and manufactured to provide uniform and flat sections every time.

**WESWOX** Advance Manual Microtome is specially designed for routine and research applications in Histopathology, Pathology, Diagnostic Labs & Medical Colleges.

# ADVANCE MANUAL ROTARY MICROTOME

*For Over 60 years,  
WESWOX MICROTOMES  
have been the classic sectioning machine  
in virtually every medical teaching and  
research laboratory and  
it's still unequalled in performance.*



## Features :

- Rust Proof Construction.
- JP Series robust disposable blade holder designed to hold either high or low profile disposable blades or both.
- Lateral knife adjustment to ensure that the entire width of the blade is used without releasing the blade clamping mechanism.
- Large storage capacity on the top of the instrument for easy access of tools.
- Finger protection guard to protect the user from any injuries caused due to blade.
- Lockable handwheel at any position allowing change of blade/specimen easier and safer.
- Robust and stable base plate for maximum stability to the microtome.
- Manual Microtome with SSC specimen feed system.
- Ergonomic and user friendly design.
- Section waste tray to collect section waste and is easily removable for cleaning purposes.

## Technical Specifications :

- Section Thickness Range : 0.5 um-60um
- Setting Values : 0.5 to 10 um in 1um increments  
10 to 20 um in 2 um increments;  
20 to 60 um in 5 um increment OR
- Minimum Setting Value : 0.5um
- Object feed : 28mm
- Vertical stroke : 60mm
- Specimen Orientation : x- and y-axes : universal 8°  
Z-axis: upto 360°
- Max. Specimen size : 50x50mm.
- Operation : Manual.
- Certification : ISO, CE approved
- Weight : 35kg.
- Dimension : 584x508x381mm

## Standard Accessories :

- High Profile Disposable blade packet of 50 no.s.
- Dust Cover.
- Section Waste Tray.
- Instruction Manual.



**D.VINODH**  
Proprietor



**VEU'S SCIENTIFIC & INSTRUMENTS**  
An Ultimate Destination For Lab Equipments  
Lab Equipment Mfg. for Research Lab,  
Engineering College, Medical College, Hospital,  
Schools etc.