



Salient Features

- ◆ Low Battery Indication
- ◆ Auto Calibration with External Weight
- ◆ Multifunction Weighing Units
- ◆ Auto Zero tracking
- ◆ Overload Protection Design
- ◆ S.S. Weighing Pan, Level indicator, Adjustable Fit
- ◆ High Resolution, Quick Weighing, Accurate Result
- ◆ In built rechargeable battery for continues use.

Precision Balance

PGB 220 & 321

HSN Code - 9016

SPECIFICATIONS	MODEL	PGB220	PGB321
	Capacity	220 g	330 g
	Readability	0.001 g (1 mg)	
	Repeatability (±)	0.001 g	
	Linearity (±)	0.002 g	
	Pan Size	90 mm Ø	
	Make	WENSAR	
	Response Time	3 - 4 Sec	
	Calibration	Automatic External	
	Display	LCD with back light	
	Tare Range	Full	
	Operating Temperature	10°C to 40°C	
	Power supply	AC 220V / 50-60 Hz	
Dimension (LxWxH) mm	360 x 260 x 355		
Weight	3.5 kg		

★ Specifications are subject to change without notice.