

Automatic External Defibrillator



Specifications

AED7000 (Automatic External Defibrillator) is a defibrillation device, which is safe, portable and easy to use. It can be applied in family and public place to give first aid for the sudden death symptom. Combined with advanced ECG analytical technique from M&B company, the defibrillator can analyze the patients ECG, make a judgement for the patients cardiac condition and indicate whether the defibrillation is needed and the defibrillation strength. It can not only enhance the success rate of defibrillation, but also reduce the injury of defibrillation to heart to the utmost extent. There is low requirement for the operation to use the device. The defibrillator will guide the operator to carry out defibrillation treatment via voice prompt step by step and give safety alert if the occasion arises. It can carry early defibrillation treatment effectively and protect the operator farthest, being very applicable to the emergency treatment for the sudden cardiac death patient in family and public place.

The AED7000 Automatic External Defibrillator delivers the defibrillation shock through two self-adhesive, pre-gelled, low-impedance electrode, defibrillator pads.

KEY FEATURES

- Three-step defibrillation process
- Two-button operation
- Extensive voice and visual prompts for the operator
- Biphasic energy output
- Lockout protection to prevent inadvertent defibrillation
- PC-based event review (optional)
- Multilingual selections are available. (optional)

