



## Ice Flaker

### LMIF SERIES

### Salient Features

- ◆ Fully Automatic microprocessor control
- ◆ Stainless Steel construction
- ◆ Continuous Ice Flakes Output
- ◆ CFC free compressor
- ◆ Low water level Indication
- ◆ Overload protection
- ◆ Noiseless operation
- ◆ Compact Design and low maintenance

### APPLICATION :

Hospitals, Educational Research Institutes & Biotechnology Laboratories

Please use Ice Flaker with stabilizer for better results & Long Life

HSN Code : 84186990

| SPECIFICATIONS       | MODEL                                  | LMIF30                          | LMIF50      | LMIF100 |
|----------------------|--|---------------------------------|-------------|---------|
|                      | Make                                   | LABMAN                          |             |         |
|                      | Capacity of Ice flake (kg per 24 hrs.) | 30 kg                           | 50 kg       | 100 kg  |
|                      | Storage Capacity (kg ~ app.)           | 5 kg                            | 7 kg        | 15 kg   |
|                      | Storage bin                            | ABS Plastic                     |             |         |
|                      | Cooling Method                         | Air                             |             |         |
|                      | Ambient temperature                    | Below 25 °C for optimum working |             |         |
|                      | Construction                           | Stainless Steel                 |             |         |
|                      | Noise level                            | ≤52 dB                          |             |         |
|                      | Water Supply                           | 3 / 4 inch BSP                  |             |         |
| Refrigerant          | R290A                                  |                                 |             |         |
| Input Power          | 180W                                   | 280W                            | 420W        |         |
| Power supply         | AC 220V / 50Hz                         |                                 |             |         |
| Dimension LxWxH (mm) | 300X493X547                            | 380X543X722                     | 548X611X883 |         |
| Gross Weight         | 36kg                                   | 45kg                            | 64kg        |         |

★ Specifications are subject to change without notice.

**STANDARD ACCESSORIES** : Main Instrument, Inlet Tube, Drainage Tube, Ice Shovel (Scoop), Instruction Manual, Test Certificate